


Anti-Cd38 Rabbit Monoclonal Antibody For Wb, If/Icc, Fc, Elisa - P56528

商品番号: CM0006482



前書き

Mouse CD38 Rabbit Monoclonal Antibody validated for WB, IF/ICC, FC, ELISA. Targets ADP-ribosyl cyclase/cyclic ADP-ribose hydrolase 1 involved in calcium signaling and lymphocyte activation.

[詳細を学ぶ](#)

Parameter	Value
Product Name	CD38 Rabbit mAb
Remarks/Alias	I-19; ADPRC 1; Cd38-rs1; CD38
Species	Mouse
Gene ID (Human)	952
Gene ID (Mouse)	12494
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 45-304 of mouse CD38 (NP_031672.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, FC, ELISA
Cross Reactivity	Mouse
SWISS	P56528
Protein Weight	34kDa
Shipping	Ice bag